

Super mega cassette*Processing consumables*

Manufacturer: Diapath S.p.A. via Savoldini 71 – 24057 Martinengo (BG) – Italy – Tel.0363-986411
 Fax 0363-948000 www.diapath.com info@diapath.com



Code	Description	Packaging
070750	White Super mega cassettes	100 pcs

Intended use

Cassettes designed for processing of big specimens, especially for brain.
 Suitable with most common microtomes on market

Composition

Acetal resin (polyacetal polyoxymethylene)

Functional features

- Rectangular grill with big holes
- Detachable double-closure system cap
- Resistant to histological reagents as alcohols, solvents, descaler
- Glazed writing area on three sides that allows specimen identification (use a permanent ink pen)
- Non-sterile
- Latex free
- Range temperature Min -40°C - Max +110°C
- Weight 4.1 g
- External size 75 x 54 x 19 mm
- Internal size 55 x 49 x 15 mm

Instructions of use

- Identify the cassette through standard procedures
- Insert the specimen into the cassette
- Close the cassette
- Proceed with standard specimen processing

To avoid mistakes, qualified and trained staff should use the product. Professional use product. The guidelines concerning safety on the workplace must be applied according to current regulations. The tools used for diagnosis must be suitable for diagnostic use in laboratory. Only authorized, trained and competent staff should perform the diagnosis

Quality control

Both products and raw materials are constantly integrated and tracked in information systems allowing their traceability.

Storage

Store the product according to the specifications listed on the label. The product, if opportunely stored and integrally packed, is stable up to the expiry date reported on the label. Do not use after expiration date.

If the product is not stored as recommended, its performance may change and must be validated by the user.

Disposal instructions

The expired and/or unused product must be disposed according to local waste regulations, based on danger classification on the label and after possible contaminations evaluation. In some cases it may be necessary an analytical evaluation to determine the correct waste classification and the danger features.

Warnings and precautions

Please refer to instructions of use.

Labeling legend



Batch n.



Manufacturer



Storage temperature limits



Product code



Expiry date



In vitro diagnostic medical device



Photosensitive