

CytoPath® Easy

Liquid-Based Cytology with no need of any other instrument



LBC It's Easy!

Easy allows a very simple cytology procedure for the liquid-based preparation of Gyn and Non-Gyn cytology specimens.

This method is intended as a replacement for the conventional Pap smear preparation for use in cervical cancer screening.

The process is simply based on gravity sedimentation and converts a liquid suspension of a cytological sample into a homogeneous thinlayer of cells. The device is composed of a specially designed vial with a cap seal press lined to be placed directly over the Cyto-Slide. Excellent results with Molecular Biology analysis on the samples stored in CytoPath® fixative. Non-toxic.

DIAPATH

DIAPATH S.p.A.

Via Savoldini,71 | 24057 Martinengo (BG) Italy | Tel. (+39) 0363.986411
Fax (+39) 0363.948000 | info@diapath.com | www.diapath.com

Follow Us



1MKTOBR160EN

CytoPath®

Complete solutions
for Liquid-Based Cytology

CytoPath®

The complete solution for an accurate, efficient and reliable cytological screening

The line consists of a complete range for cyto specimens:



Medscand Pap Perfect Spatula



Medscand Cytobrush

+



Rovers Cervex Brush

COLLECTION
DEVICES



CytoPath® Gyn



CytoPath® NON-Gyn



CytoPath® FISH

CONSUMABLES



For high quality performances

CONTAINERS



For brilliant and consistent results

STAINING
LINE

CytoPath® Processor

Significantly more effective

CytoPath® Processor is a new instrument processor specific for the preparation of liquid-based cytological specimen



Confidence in results

The semi-automatic instrument allows processing both GYN and NON-GYN samples. All available consumables are in CytoPath® line and consist of cytological fixative, filters (GYN and NON-GYN according to samples) and slides.

Increase efficiency

Processing samples using CytoPath® Processor allows obtaining slide with thinlayer cytological preparation in spots of 20 mm of diameter.

